








Histamine and 1-methylhistamine in Serum LC-MS/MS Analysis Kit

Histamine intolerance is a disorder of histamine homeostasis due to decreased intestinal degradation of this amine; it is mainly caused by a deficiency of the enzyme diamine oxidase, which leads to its accumulation in plasma and the development of adverse health effects. Two metabolic pathways are known for histamine in humans. Histamine-N-methyltransferase is the enzyme responsible for the ring methylation of histamine, found primarily in the liver and kidney, and ultimately converts histamine to 1-methylhistamine.

Highlights of the Analysis Kit

-  Just a few pipetting steps for the sample preparation; for urine “**protein crash and shoot**”
-  Total run time is 6.5 min.
-  Use of stable isotope labeled internal standard
-  Small volume of patient’s sample is required
-  Long life span of HPLC column

Parameters

Histamine and 1-methylhistamine


Sample Type

Serum

Sample Preparation

| | |
|---|--|
| 1 | Pipette 50 μ L of calibrator/control/serum sample into a glass centrifuge tube |
| 2 | Then, add 50 μ L of internal standard and 600 μ L of Reagent-1 respectively, vortex 5 sec. then, centrifuge at 4000 rpm for 5 min. |
| 3 | Decant the clear supernatant into a HPLC vial prior to injection |

Example Chromatogram



Total ion chromatogram of histamine and 1-methylhistamine

Method Performance

All validation results were obtained using Agilent 6465 TQ (Ultivo) system

| Analyte | LOQ (ng/ml) | Linearity (ng/ml) | Recovery | | Repeatability | | | |
|-------------------|----------------|----------------------|--------------|---------------|---------------|---------------|---------------|---------------|
| | | | LLQC* (%) | HLQC** (%) | intra-day | | inter-day | |
| | | | | | LLQC (%CV) | HLQC (%CV) | LLQC (%CV) | HLQC (%CV) |
| Histamine | 0.26 | 0.5 – 25.0 | 104 | 103 | 5.41 | 2.07 | 5.74 | 3.40 |
| 1-Methylhistamine | 0.23 | 0.5 – 25.0 | 102 | 111 | 4.68 | 2.64 | 4.82 | 3.69 |

* LLQC: Low-level quality control

** HLQC: High-level quality control



Altium International Laboratuvar Cihazları A.Ş.
Barbaros Mah. Temmuz Sk. No:6 Altium Plaza Ataşehir, İstanbul
T: +90 216 571 02 00 F: +90 216 571 02 02

www.jasem.com.tr